

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human NPC2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12129
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human NPC2 (Asp31-Leu151).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HE1, NPC intracellular cholesterol transporter 2, He1, Epididymal secretory protein E1, Human epididymis-specific protein 1, NPC2, Niemann-Pick disease type C2 protein