

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TCF12 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12137
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human TCF12 (Asn567-Ser658).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-bHLHb20, HTF4, E-box-binding protein, HEB, TCF12, Class B basic helix-loop-helix protein 20, Transcription factor 12, BHLHB20, DNA-binding protein HTF4, Transcription factor HTF-4, TCF-12