

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human RBCK1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12165
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human RBCK1 (Met1-Gly156).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-RING finger protein 54, RNF54, RanBP-type and C3HC4-type zinc finger-containing protein 1, C20orf18, XAP4, Hepatitis B virus X-associated protein 4, UBCE7IP3, RBCK1, RING-type E3 ubiquitin transferase HOIL-1, HOIL-1, Heme-oxidized IRP2 ubiquitin liga