

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human EIF3B Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12181
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human EIF3B (Arg303-Arg707).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-eIF3 p110, Eukaryotic translation initiation factor 3 subunit B, eIF-3-eta, eIF3b, EIF3B, Prt1 homolog, Eukaryotic translation initiation factor 3 subunit 9, EIF3S9, hPrt1, eIF3 p116