

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human GPAT3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12198
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human GPAT3 (Ser207-Ser434).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-1-acyl-sn-glycerol-3-phosphate O-acyltransferase 9, AGPAT 10, 1-AGPAT 9, GPAT-3, hGPAT3, Glycerol-3-phosphate acyltransferase 3, Lysophosphatidic acid acyltransferase theta, 1-acyl-sn-glycerol-3-phosphate O-acyltransferase 10, LPAAT-theta, Acyl-CoA:g