

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human MAP1LC3C Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12221
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human MAP1LC3C (Arg11-Leu147).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Microtubule-associated proteins 1A/1B light chain 3C, MAP1A/MAP1B LC3 C, Autophagy-related ubiquitin-like modifier LC3 C, MAP1A/MAP1B light chain 3 C, Microtubule-associated protein 1 light chain 3 gamma, MAP1 light chain 3-like protein 3, MAP1LC3C,