

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TRIM15 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12224
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human TRIM15 (Glu80-Pro275).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ZNF178, Tripartite motif-containing protein 15, TRIM15, ZNFB7, Zinc finger protein 178, RNF93, RING finger protein 93, Zinc finger protein B7