

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CBX1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12291
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human CBX1 (Glu24-Leu168).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-HP1Hsbeta, CBX, Chromobox protein homolog 1, Modifier 1 protein, Heterochromatin protein 1 homolog beta, HP1 beta, M31, CBX1, p25beta, Heterochromatin protein p25