

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human AP2A1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12334
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human AP2A1 (Ile737-Phe977).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-AP-2 complex subunit alpha-1, Alpha1-adaptin, Adaptor protein complex AP-2 subunit alpha-1, 100 kDa coated vesicle protein A, Clathrin assembly protein complex 2 alpha-A large chain, ADTAA, CLAPA1, AP2A1, Adaptor-related protein complex 2 subunit alp