

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human U2AF1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12339
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human U2AF1 (Met1-Arg196).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-U2AFBP, U2 auxiliary factor 35 kDa subunit, U2AF1, U2 snRNP auxiliary factor small subunit, U2 small nuclear RNA auxiliary factor 1, Splicing factor U2AF 35 kDa subunit, U2AF35