

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human NFIL3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12354
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human NFIL3 (Glu65-Ser161).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-IL3BP1, Nuclear factor interleukin-3-regulated protein, E4 promoter-binding protein 4, Interleukin-3 promoter transcriptional activator, NFIL3, Interleukin-3-binding protein 1, Transcriptional activator NF-IL3A, E4BP4