

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human RPSA Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12362
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human RPSA (Met1-Arg186).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-LAMBR, 37LRP, LBP/p40, Multidrug resistance-associated protein MGr1-Ag, Colon carcinoma laminin-binding protein, LamR, 67 kDa laminin receptor, NEM/1CHD4, Laminin receptor 1, Small ribosomal subunit protein uS2, Laminin-binding protein precursor p40,