

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Porphyrromonas gingivalis kgp/prtK Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12386
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Porphyromonas gingivalis kgp/prtK (Asp229-Pro683).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Lys-gingipain W83, 3.4.22.47, Lysine specific cysteine protease, Lysine-specific cysteine proteinase, Porhypain, PrtK48, Lys-gingipain catalytic subunit, 39 kDa adhesin, PrtK39, 15 kDa adhesin, PrtK15, 44 kDa adhesin, PrtK44, kgp, prtK, prtP