

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ATF6B Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12403
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human ATF6B (Ala2-Arg50).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-CREBL1, Protein G13, cAMP-responsive element-binding protein-like 1, Activating transcription factor 6 beta, cAMP response element-binding protein-related protein, Cyclic AMP-dependent transcription factor ATF-6 beta, G13, ATF6B, Creb-rp, ATF6-beta,