

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human UHRF1 (Gly427-Gly636) Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12472
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human UHRF1 (Gly427-Gly636).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Nuclear zinc finger protein Np95, RING-type E3 ubiquitin transferase UHRF1, RING finger protein 106, Transcription factor ICBP90, E3 ubiquitin-protein ligase UHRF1, Ubiquitin-like PHD and RING finger domain-containing protein 1, NP95, HuNp95, Nuclear