

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human RUVBL2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12489
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human RUVBL2 (Gln49-Ser463).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-48 kDa TATA box-binding protein-interacting protein, RuvB-like 2, INO80 complex subunit J, RUVBL2, TAP54-beta, Reptin 52, TIP60-associated protein 54-beta, TIP48, TIP49B, 51 kDa erythrocyte cytosolic protein, 48 kDa TBP-interacting protein, Repressin