

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human PMEL Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12508
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human PMEL (Gln65-Thr299).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ME20M, Melanocyte protein PMEL, Silver locus protein homolog, Secreted melanoma-associated ME20 antigen, D12S53E, 95 kDa melanocyte-specific secreted glycoprotein, SILV, P100, Premelanosome protein, PMEL17, ME20S, PMEL, ME20-S, ME20-M, P26, P1, Melan