

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human HSPA9/Mortalin/GRP75 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12510
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human HSPA9/Mortalin/GRP75 (Lys52-Asp431).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-mt-HSP70, HSPA9B, MOT, Peptide-binding protein 74, Stress-70 protein, mitochondrial, PBP74, GRP-75, Mortalin, Heat shock 70 kDa protein 9, GRP75, 75 kDa glucose-regulated protein, HSPA9