

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human ZC3HAV1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12541
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human ZC3HAV1 (Met1-Asn227).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Zinc finger CCCH-type antiviral protein 1, Zinc finger CCCH domain-containing protein 2, ZC3HAV1, ZC3HDC2, ZAP, Zinc finger antiviral protein, ARTD13, PARP13, ADP-ribosyltransferase diphtheria toxin-like 13, Inactive Poly [ADP-ribose] polymerase 13