

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human UQCRFS1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12560
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human UQCRFS1 (Arg104-Gly274).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Rieske iron-sulfur protein, UQCRFS1 mitochondrial targeting sequence, UQCRFS1, Ubiquinol-cytochrome c reductase 8 kDa protein, Cytochrome b-c1 complex subunit 11, Subunit 9, Su9, RISP, Rieske protein UQCRFS1, Cytochrome b-c1 complex subunit 5, Comple