

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TRIM63 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12621
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human TRIM63 (Gly11-Arg112).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-IRF, RNF28, Iris RING finger protein, SMRZ, MURF1, MuRF-1, MuRF1, Striated muscle RING zinc finger protein, Tripartite motif-containing protein 63, TRIM63, E3 ubiquitin-protein ligase TRIM63, Muscle-specific RING finger protein 1, RING finger protein