

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CAPRIN1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12628
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human CAPRIN1 (Arg132-Asn251).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-GPI-anchored protein p137, Cytoplasmic activation- and proliferation-associated protein 1, p137GPI, RNG105, M11S1, Caprin-1, Cell cycle-associated protein 1, RNA granule protein 105, CAPRIN1, GPIP137, GPI-anchored membrane protein 1, GPI-p137, GPIAP1