

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human OPTN Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12650
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human OPTN (Ser415-Arg524).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-NRP, TFIIIA-IntP, Huntingtin-interacting protein 7, HIP-7, Huntingtin yeast partner L, Transcription factor IIIA-interacting protein, NEMO-related protein, E3-14.7K-interacting protein, GLC1E, HIP7, HYPL, FIP-2, Huntingtin-interacting protein L, Opti