

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Rat SCNN1A/SCNN1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12679
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Rat SCNN1A/SCNN1 (Glu132-Asp391).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Rattus norvegicus (Rat)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Amiloride-sensitive sodium channel subunit alpha, Alpha-NaCH, Epithelial Na(+) channel subunit alpha, Alpha-ENaC, Nonvoltage-gated sodium channel 1 subunit alpha, SCNEA, Scnn1a, Renac