

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human RACK1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12742
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human RACK1 (Met1-Arg317).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-RACK1, Human lung cancer oncogene 7 protein, Guanine nucleotide-binding protein subunit beta-like protein 12.3, GNB2L1, Cell proliferation-inducing gene 21 protein, Guanine nucleotide-binding protein subunit beta-2-like 1, Receptor of activated prote