

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TNIP1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12777
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human TNIP1 (Glu2-Gln636).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Naf1, Virion-associated nuclear shuttling protein, VAN, TNFAIP3-interacting protein 1, KIAA0113, HIV-1 Nef-interacting protein, hVAN, A20-binding inhibitor of NF-kappa-B activation 1, Nef-associated factor 1, TNIP1, NAF1, Nip40-1, ABIN-1