

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-DENV-3 Protein prM/prM Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12782
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant DENV-3 Protein prM/prM (Gly122-Arg205).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Genome polyprotein, Capsid protein C, Core protein, Protein prM, Peptide pr, Small envelope protein M, Matrix protein, Envelope protein E, Non-structural protein 1, NS1, Non-structural protein 2A, NS2A, Serine protease subunit NS2B, Flavivirin protea