

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse IL36RN/IL-1F5/IL-36Ra Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12791
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Mouse IL36RN/IL-1F5/IL-36Ra (Met1-Asp156).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-FIL1D, Interleukin-1 family member 5, IL-36Ra, IL1HY1, Interleukin-36 receptor antagonist protein, IL1RP3, IL-1RP3, IL-1ra homolog 1, Interleukin-1-like protein 1, IL-1-related protein 3, IL1L1, Interleukin-1 HY1, IL-1HY1, Interleukin-1 delta, FIL1 d