

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse GIT2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12795
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Mouse GIT2 (Phe29-Gly358).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-KIAA0148, Cool-interacting tyrosine-phosphorylated protein 2, CAT-2, GRK-interacting protein 2, GIT2, G protein-coupled receptor kinase-interactor 2, ARF GTPase-activating protein GIT2, ARF GAP GIT2, CAT2