

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human PHLDA1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A12955
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	E. coli - derived recombinant Human PHLDA1 (Ser137-Gln305).

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Proline- and histidine-rich protein, T-cell death-associated gene 51 protein, Pleckstrin homology-like domain family A member 1, PQR protein, Apoptosis-associated nuclear protein, PHRIP, PHLDA1, Proline- and glutamine-rich protein, TDAG51, PQ-rich pr