

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Mouse VISTA & PD-L1 Bispecific Antibody (lv0246)
PRODUCT REFERENCE	ARO-A12974
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Bispecific

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bispecific, V-type immunoglobulin domain-containing suppressor of T-cell activation, Platelet receptor Gi24, V-set domain-containing immunoregulatory receptor, V-set immunoregulatory receptor, Vsir, Dies1, PD-1H, VISTA, Programmed cell death 1 ligand