

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Mouse PD-1 & CD47 Bispecific Antibody (Iv0241)
PRODUCT REFERENCE	ARO-A12979
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Bispecific

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bispecific, Programmed cell death protein 1, Protein PD-1, mPD-1, CD279, Pdcd1, Pd1, Leukocyte surface antigen CD47, Integrin-associated protein, IAP, CD47, Cd47