

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	InVivoMAb Anti-hERBB2/HER-2 & mCD3e Bispecific Antibody (Iv0230)
<b>PRODUCT REFERENCE</b>	ARO-A12987
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	Bispecific

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Bispecific, Receptor tyrosine-protein kinase erbB-2, 2.7.10.1, Metastatic lymph node gene 19 protein, MLN 19, Proto-oncogene Neu, Proto-oncogene c-ErbB-2, Tyrosine kinase-type cell surface receptor HER2, p185erbB2, CD340, ERBB2, HER2, MLN19, NEU, NGL