

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-hEGFRvIII & mCD3e Bispecific Antibody (Iv0229)
PRODUCT REFERENCE	ARO-A12988
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Bispecific

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bispecific, Epidermal growth factor receptor, 2.7.10.1, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1, EGFR, ERBB, ERBB1, HER1, T-cell surface glycoprotein CD3 epsilon chain, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, Cd3