

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Mouse CTLA4 & OX40 Bispecific Antibody (lv0227)
PRODUCT REFERENCE	ARO-A12990
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Bispecific

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bispecific, Cytotoxic T-lymphocyte protein 4, Cytotoxic T-lymphocyte-associated antigen 4, CTLA-4, CD152, Ctla4, Cd152, Tumor necrosis factor receptor superfamily member 4, OX40 antigen, OX40L receptor, CD134, Tnfrsf4, Ox40, Txgp1