

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Streptococcus pneumoniae iga/IgA1 protease Antibody (mAb2)
PRODUCT REFERENCE	ARO-A12999
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Immunoglobulin A1 protease

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Immunoglobulin A1 protease, IgA1 protease, 3.4.24.13, IgA-specific zinc metalloproteinase, iga, spr1042