

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME InVivoMAb Anti-Human BMP4 Antibody (6H4)

PRODUCT REFERENCE ARO-A13011

CLONALITY Monoclonal Antibody

ISOTYPE

TARGET Bone morphogenetic protein 4

PRODUCT INFORMATION

ORIGIN SPECIES *Mus musculus (Mouse)*

FORM Liquid

BUFFER 0.01M PBS (pH 7.4)

STABILITY & STORAGE Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS Anti-Bone morphogenetic protein 4, BMP-4, Bone morphogenetic protein 2B, BMP-2B, BMP4, BMP2B, DVR4