

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-DENV-1/2/3/4 Envelope protein E/DII-FL region Antibody (1C19)
PRODUCT REFERENCE	ARO-A13038
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Genome polyprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Genome polyprotein, Envelope protein E, DII-FL region, Dengue virus type 1 (strain Nauru/West Pac/1974) (DENV-1), Dengue virus type 2 (strain Jamaica/1409/1983) (DENV-2), Dengue virus type 3 (strain Philippines/H87/1956) (DENV-3), Dengue virus type 4