

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-DENV-2 Envelope protein E/EDIII domain Antibody (3H5)
PRODUCT REFERENCE	ARO-A13039
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Genome polyprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Genome polyprotein, Envelope protein E, EDIII domain, Dengue virus 2, Dengue virus type 2 (strain Thailand/PUO-218/1980) (DENV-2)