

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME InVivoMAb Anti-LASV GPC Antibody (LAVA01)

PRODUCT REFERENCE ARO-A13049

CLONALITY Monoclonal Antibody

ISOTYPE

TARGET Pre-glycoprotein polyprotein GP complex

PRODUCT INFORMATION

ORIGIN SPECIES *Mus musculus (Mouse)*

FORM Liquid

BUFFER 0.01M PBS (pH 7.4)

STABILITY & STORAGE Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS Anti-Pre-glycoprotein polyprotein GP complex, Pre-GP-C, Stable signal peptide, SSP, Glycoprotein G1, GP1, Glycoprotein G2, GP2, GPC, GP-C, Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV)