

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	InVivoMAb Anti-DENV2 Protein prM/prM Antibody (2H2)
<b>PRODUCT REFERENCE</b>	ARO-A13090
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	Genome polyprotein

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-Genome polyprotein, Protein prM, Peptide pr, Small envelope protein M, Matrix protein, Envelope protein E, Non-structural protein 1, NS1, Dengue virus type 2 (strain Thailand/PUO-218/1980) (DENV-2), Dengue virus type 2 (strain Jamaica/1409/1983) (DENV)