

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-ZEBOV GP/Envelope glycoprotein Antibody (6D8)
PRODUCT REFERENCE	ARO-A13102
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Envelope glycoprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Envelope glycoprotein, GP1,2, GP, Shed GP, GP1,2-delta, GP1, GP2, GP, Zaire ebolavirus (strain Mayinga-76) (ZEBOV) (Zaire Ebola virus)