

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Escherichia coli shiga toxin type 2/stx2A Antibody (11E10#)
PRODUCT REFERENCE	ARO-A13106
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	rRNA N-glycosylase

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-rRNA N-glycosylase, 3.2.2.22, stx2 A-subunit, stx-2 A subunit, stx2, stx2A, stx2a, stx2A_1, stx2cA, stx2vA, stxA2, stxA_1, stxA_3, vta, vtx2A, Escherichia coli