

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-Escherichia coli shiga toxin type 2/stx2B Antibody (lv0194)
PRODUCT REFERENCE	ARO-A13107
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Shiga toxin 2 B subunit

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Shiga toxin 2 B subunit, Shiga toxin 2B, Shiga toxin 2c subunit B, Shiga toxin 2d subunit B, Shiga toxin 2v subunit B, Shiga-like toxin II B-subunit, Shiga-like toxin type IIvhc, Stx2c protein B subunit, VT2va beta-subunit, Verocytotoxin 2d B-subunit