

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	InVivoMAb Anti-Escherichia coli shiga toxin type 1/stx1B Antibody (lv0193)
<b>PRODUCT REFERENCE</b>	ARO-A13108
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	SLT-I B Subunit

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-SLT-I B Subunit, Shiga toxin 1 B-subunit, Shiga toxin 1c subunit B, Shiga toxin Stx1 subunit B, Shiga toxin Stx1c subunit B, Shiga-like toxin 1B-subunit, Sxt1B, stx1vB, stx1 B-subunit, stx1B, stx1b, stxB1c, Escherichia coli