

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | InVivoMAb Anti-HBV-A S/L-HBsAg/L glycoprotein Antibody (Iv0190) |
| PRODUCT REFERENCE | ARO-A13112 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | |
| TARGET | Large envelope protein |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Mus musculus (Mouse)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-Large envelope protein, L glycoprotein, L-HBsAg, LHB, Large S protein, Large surface protein, Major surface antigen, S, Hepatitis B virus genotype A2 subtype adw (isolate Japan/Nishioka/1983) (HBV-A) |