

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-HBV-A S/L-HBsAg/L glycoprotein Antibody (Iv0189)
PRODUCT REFERENCE	ARO-A13113
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Large envelope protein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Large envelope protein, L glycoprotein, L-HBsAg, LHB, Large S protein, Large surface protein, Major surface antigen, S, Hepatitis B virus genotype A2 subtype adw (isolate Japan/Nishioka/1983) (HBV-A)