

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-HHV-5/HCMV Protein UL128 & UL131A Antibody (lv0187)
PRODUCT REFERENCE	ARO-A13119
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Uncharacterized protein UL128

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Uncharacterized protein UL128, UL128, Protein UL131A, UL131A, Human cytomegalovirus (strain AD169) (HHV-5) (Human herpesvirus 5), Human cytomegalovirus (strain Merlin) (HHV-5) (Human herpesvirus 5)