



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME InVivoMAb Anti-Human TGM2 Antibody (52F3)

PRODUCT REFERENCE ARO-A13121

CLONALITY Monoclonal Antibody

ISOTYPE

TARGET Protein-glutamine deamidase TGM2

PRODUCT INFORMATION

ORIGIN SPECIES *Mus musculus (Mouse)*

FORM Liquid

BUFFER 0.01M PBS (pH 7.4)

STABILITY & STORAGE Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS Anti-Protein-glutamine deamidase TGM2, TGase II, Heart G alpha(h), Erythrocyte transglutaminase, Protein-glutamine gamma-glutamyltransferase 2, TG(C), Transglutaminase C, Protein-glutamine serotonyltransferase TGM2, tTG, tTgase, Transglutaminase-2, Protei

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY