

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	InVivoMAb Anti-ZIKV/DENV NS1c Antibody (1G5)
PRODUCT REFERENCE	ARO-A13127
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Genome polyprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Genome polyprotein, Capsid protein C, Capsid protein, Core protein, Non-structural protein 1, NS1, C-terminal fragment of NS1, NS1c, Zika virus (isolate ZIKV/Human/French Polynesia/10087PF/2013) (ZIKV), Dengue virus type 2